TRACOE[®] purofoam Tracheal Compress

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Instructions for Use for TRACOE purofoam Tracheal Compress

Note: Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your own and your patients' safety, please follow instructions below.

1. Intended Use and Indications for Use

The TRACOE purofoam Tracheal Compress is designed to help protecting the skin surrounding and wound environment of a tracheostoma. Its main functions are the following:

- absorbing secretions from the tracheostoma
- acting as a mechanical barrier between the skin and the neck flange of the tracheostomy tube
- providing a waterproof but breathable (vapor permeable) barrier, that also aids in the prevention of external bacterial contamination

Patient Population: neonates, infants, children, adolescents and adults

Intended User: The product is intended for mechanically ventilated and self-breathing patients in hospitals, pre-hospitals (EMS), extended care facilities, outpatient clinics and can be used by trained individuals in tracheostomy care.

Indications: The TRACOE purofoam Tracheal Compress is used for patients wearing a tracheostomy tube.


Single Patient Use and Useful Life: This is a single use disposable product. The change intervals can reach up to 3 days, depending on the amount of exudation or secretions. If the compress is very swollen or leakage occurs the compress should be changed more frequently.

2. General Description

The TRACOE purofoam Tracheal Compress is made of an absorbent open pore polyurethane foam that maintains the appropriate moisture to wound environment. The surface layer of the compress that faces the skin allows the absorption of liquids into the foam's core. This side is covered with a peelable protective foil until use. The TRACOE purofoam Tracheal Compress can absorb secretions and exudate coming from the stoma. The outside layer is vapor permeable and liquid-proof. It also acts as an additional mechanical,

microbiological barrier. The compress meets modern standards of moist wound treatment.

Due to the circular opening and the zigzag slit, the compress can be placed or removed behind the neck flange of the tracheostomy tube while it is in place. A crescent shaped recess on one side of the compress allows placement under the chin. The low profile of the compress helps to make the compress comfortable for the patient.

The TRACOE purofoam Tracheal Compress is provided individually packaged in a sterile pouch.

3. Contraindications

- Do not use the compress for festering wounds.
- The product must not be used in case of known sensitivity to plasters, adhesive tapes, or wetting agents.

4. General Precautions

- When the product is used together with other medical devices, follow their respective instructions for use. Contact the manufacturer if there are any questions, or if assistance is required.
- Safety precautions must be taken in the event of complications during the described procedures.
- The product should be inspected for integrity and function prior to use. If the product is damaged, it should be replaced with a new product.
- The use of ointments, balms and wetting agents may damage or destroy parts of the compress as well as cleaning / disinfection with hypochlorite solution or oxygenated water.

5. Warnings

- Do not use this product if the sterile packaging has been compromised or damaged, as this may impair the sterility and integrity of the product.
- Refurbishing (including re-sterilisation) is not allowed, this may influence the material and function of the product. The product is single use only.

6. Adverse Reactions

- Longer application times of the compress may lead to a rash.
- If allergic reactions such as itching, redness, or other signs of hypersensitivity occur, remove the compress immediately and consult a physician.

7. Functional Description

7.1 Preparation

The following functions must be checked immediately prior to use. If the device fails the initial inspection, repeat the procedure with a new device. Do not discard the device, and follow instructions provided in section "Returns and Complaints".

1. Inspect the sterile packaging to ensure it is secure, undamaged and all components are present.
2. Open the package and inspect the device for damage prior to use.

7.2 Preparing the Patient

1. Clean and dry the tracheostoma and the surrounding skin.
2. Skin and tracheostoma should be free of disinfectants or cleaning agents.

7.3 Procedure

Caution:

- Do not dislocate the tracheostomy tube when applying or removing the compress.
- Do not disconnect the gas supply during the application.
- TRACOE purofoam Tracheal Compress should be inspected frequently at short intervals.

1. Choose the best fitting compress size depending on the tube diameter and the amount of secretion/exudate.
2. If necessary, the compress may be cut to a specific shape.
3. Use the compress immediately after opening the sterile pouch.
4. Remove the protective foil on the patient side of the compress.
5. Gently lift the neck flange of the tracheostomy tube to create a small gap for the compress.
6. Place the compress around the tube using the zigzag slit.
7. The orientation of the compress can be adapted to the local conditions by slight rotation.
8. The crescent shaped recess should be below the chin.
9. If necessary, the compress can be fixed with additional plasters or medical adhesive tapes. However, avoid covering the compress completely. This will affect its permeability.

7.4 Changing the Compress

1. To remove the compress, lift the neck flange slightly and carefully pull the compress away from the tube.
2. If the compress does not loosen, apply standard wound treatment measures for compress exchange.
3. Position a new compress as described in chapter "Procedure".

8. Care and Cleaning

The product is single-use and disposable, therefore, no cleaning is required or allowed.

9. Storage

Store the TRACOE products in their original packaging according to the conditions displayed on the packaging.

10. Packaging

The product is provided sterile (by irradiation) which allows application under sterile conditions.

11. Disposal

Used products are to be disposed of in accordance with national regulations, waste management plans, or clinical procedures governing biohazardous waste materials, e.g. the direct disposal in a tear and moisture-resistant and secure bag or container, which is routed to the local waste disposal system for contaminated medical products.

For further recommendations, contact your hygiene officer in health facilities, or the local waste management for homecare use.

12. Returns and Complaints

Returned products, that have been used, will only be accepted, if TRACOE has agreed to the return and a completed decontamination certificate and complaint report is enclosed. These forms are available either directly from TRACOE medical, or via the website www.tracoe.com.

If the device is involved in a serious incident, as defined in Medical Device Regulation (EU) 2017/745, you may contact TRACOE medical (complaints@tracoe.com), and the Competent Authority in the country of use.

13. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be affected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either from TRACOE medical GmbH or on our website at www.tracoe.com.

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